



**GUNZE**

ABSORBABLE POLYGLYCOLIC ACID FELT


**NEOVEIL™**



# **NEOVEIL™** Staple Line Reinforcement


**FOR OPEN SURGERY & ENDOSCOPY**

**NEOVEIL™** Sheet




## **What's NEOVEIL™?**

- 100% Polyglycolic Acid
- Bioabsorbable
- Reinforcement of Soft Tissues / Organs



## Usable all over the human body



ENT, Oral & Maxillofacial surgery

Thoracic surgery

Hepatobiliary & Pancreatic surgery

ESD (Endoscopic submucosal dissection)

### Indications For Use


*Indicated for use in surgical procedures in which soft tissue transection or resection with staple or suture line reinforcement is needed. NEOVEIL™ can be used for reinforcement of suture or staple lines during lung resection, liver resection, bronchial, bariatric, colon, colorectal, esophagus, gastric, mesentery, pancreas, and small bowel procedures.*

NEOVEIL™ has been used successfully in various areas of the human body. Ask the GUNZE Team for case reports to get more information. Make sure to consult the IFU before application.

## Retention of strength

### NEOVEIL™ degrades by hydrolysis

- Completely absorbed after about 15 weeks



## Did you know...

**NORM**

...that NEOVEIL™ Sheet is the norm in Japan for preventing air leakages?

**2.61 MILLION**

Total sales from 1992 to Oct. 2021 = **2,610,000 pieces**

**0.0003%**

Rate of complaints reported to the authorities was **0.000268%** = **7 cases from 1992** for NEOVEIL™ SLR and 0 cases for NEOVEIL™ Sheet


\*According to GUNZE PMS data



## PRODUCT VARIATION, NEOVEIL™ Sheet



| Product Code     | Size         |
|------------------|--------------|
| <b>NV-L-015G</b> | 100 x 100 mm |
| <b>NV-M-015G</b> | 100 x 50 mm  |
| <b>NV-L-03G</b>  | 100 x 100 mm |
| <b>NV-M-03G</b>  | 100 x 50 mm  |
| <b>NV-L-05G</b>  | 100 x 100 mm |
| <b>NV-M-05G</b>  | 100 x 50 mm  |

For more information, please contact at [inquiry\\_gius@gunze.co.jp](mailto:inquiry_gius@gunze.co.jp)  
Products listed may not be available in all markets.



1 sheet / bag

### Full scale




## PRODUCT VARIATION, NEOVEIL™ Staple Line Reinforcement

Product Code

**NV-ET-M60E-3**

**NV-ET-M60L-3**

**NV-ET-M60A-2**



1 set ( 2 sleeves) / package

### Product Code & Compatible stapler models

| Product Code        | Compatible stapler models   |
|---------------------|---|
| <b>NV-ET-M60E-3</b> | [Johnson & Johnson]<br>ECHELON ENDOPATH 60 mm Stapler<br>ECHELON FLEX ENDOPATH 60 mm Stapler<br>ECHELON FLEX Powered ENDOPATH 60 mm Stapler   |
|                     | [Medtronic]<br>ENDO GIA II 60 Cartridge (Green)<br>ENDO GIA Reticulator 60 Cartridge (Green)<br>ENDO GIA 60 mm Articulating Extra-Thick Reload with Tri-Staple Technology (Black)   |
|                     | [Lexington Medical]<br>AEON Endoscopic Stapler Reloads AESR60B  |
| <b>NV-ET-M60L-3</b> | [Lexington Medical]<br>AEON Endoscopic Stapler Reloads AESR60P  |
| <b>NV-ET-M60A-2</b> | [Medtronic]<br>ENDO GIA II 60 Cartridge (Blue, White)<br>ENDO GIA Reticulator 60 Cartridge (Blue, White)<br>ENDO GIA 60 mm Articulating Medium/Thick Reload with Tri-Staple Technology (Purple)<br>ENDO GIA 60 mm Articulating Vascular/Medium Reload with Tri-Staple Technology (Tan)<br>ENDO GIA 60 mm Articulating Vascular Reload with Tri-Staple Technology (Gray) |

The color in parentheses is the color of the cartridge attached to the stapler. It is not the color of NEOVEIL™.

**NEOVEIL™ Sheet**  
Product video



**NEOVEIL™ SLR**  
Product video



Distributor

# GUNZE

**GUNZE LIMITED Medical Division**

<http://www.gunze.co.jp/e/medical/>

46, Natsumegaichi, Aono-cho, Ayabe,  
Kyoto 623-8513, Japan

**U.S.A.**

**GUNZE INTERNATIONAL USA, INC.**

390 Fifth Avenue, Suite 903, New York,  
NY 10018, U.S.A.

Tel: +1-212-354-9060

Fax: +1-212-354-6171

[inquiry\\_gius@gunze.co.jp](mailto:inquiry_gius@gunze.co.jp)

<https://gunzemedicalusa.com/>

**Europe**

**GUNZE INTERNATIONAL  
EUROPE GmbH**

Oststraße 80, 40210 Düsseldorf, Germany

Tel: +49-211-3613-971

Fax: +49-211-3613-164

**China**

**GUNZE MEDICAL DEVICES  
(SHENZHEN) LIMITED**

Room 802, Kerry Centre, 2008 Renminnan Rd.,  
SHENZHEN 518001 P.R. China

Tel: +86-755-8230-0553

Fax: +86-755-8230-0469

**Japan/Asia Pacific**

**GUNZE MEDICAL LIMITED  
TOKYO OFFICE**

15F, Tokyo Shiodome Bldg., 1-9-1,

Higashi Shimbashi, Minato-ku,

Tokyo 105-7315, Japan

Tel: +81-3-4485-0007

Fax: +81-3-4485-0040