

ABSORBABLE POLYGLYCOLIC ACID FELT

NEOVEIL™






NEOVEIL™

Staple Line Reinforcement


FOR OPEN SURGERY & ENDOSCOPY

NEOVEIL™ Sheet




What's NEOVEIL™?

- 100% Polyglycolic Acid
- Bioabsorbable
- Reinforcement of Soft Tissues / Organs



Usable all over the human body



ENT, Oral & Maxillofacial surgery

Thoracic surgery

Hepatobiliary & Pancreatic surgery

ESD (Endoscopic submucosal dissection)

Indications For Use


Indicated for use in surgical procedures in which soft tissue transection or resection with staple or suture line reinforcement is needed. NEOVEIL™ can be used for reinforcement of suture or staple lines during lung resection, liver resection, bronchial, bariatric, colon, colorectal, esophagus, gastric, mesentery, pancreas, and small bowel procedures.

NEOVEIL™ has been used successfully in various areas of the human body. Ask the GUNZE Team for case reports to get more information. Make sure to consult the IFU before application.

Retention of strength

NEOVEIL™ degrades by hydrolysis

- Completely absorbed after about 15 weeks



Did you know...

NORM

...that
NEOVEIL™ Sheet is the norm in Japan for preventing air leakages?

**2.61
MILLION**

Total sales from
1992 to Oct. 2021
**=2,610,000
pieces**

0.0003%


Rate of complaints
reported to the
authorities was
0.000268%
=7 cases from 1992
for NEOVEIL™ SLR and
0 cases for NEOVEIL™ Sheet

*According to GUNZE PMS data

PRODUCT VARIATION, NEOVEIL™ Sheet



Product Code	Size
NV-L-015G	100 x 100 mm
NV-M-015G	100 x 50 mm
NV-L-03G	100 x 100 mm
NV-M-03G	100 x 50 mm

For more information, please contact at inquiry_gius@gunze.co.jp
Products listed may not be available in all markets.



1 sheet / bag

Full scale




PRODUCT VARIATION, NEOVEIL™ Staple Line Reinforcement

Product Code

NV-ET-M60E-3

NV-ET-M60L-3

NV-ET-M60A-2



1 set (2 sleeves) / package

Product Code & Compatible stapler models

Product Code	Compatible stapler models
NV-ET-M60E-3	[Johnson & Johnson] ECHELON ENDOPATH 60 mm Stapler ECHELON FLEX ENDOPATH 60 mm Stapler ECHELON FLEX Powered ENDOPATH 60 mm Stapler
	[Medtronic] ENDO GIA II 60 Cartridge (Green) ENDO GIA Reticulator 60 Cartridge (Green) ENDO GIA 60 mm Articulating Extra-Thick Reload with Tri-Staple Technology (Black)
	[Lexington Medical] AEON Endoscopic Stapler Reloads AESR60B
	[Lexington Medical] AEON Endoscopic Stapler Reloads AESR60P
NV-ET-M60A-2	[Medtronic] ENDO GIA II 60 Cartridge (Blue, White) ENDO GIA Reticulator 60 Cartridge (Blue, White) ENDO GIA 60 mm Articulating Medium/Thick Reload with Tri-Staple Technology (Purple) ENDO GIA 60 mm Articulating Vascular/Medium Reload with Tri-Staple Technology (Tan) ENDO GIA 60 mm Articulating Vascular Reload with Tri-Staple Technology (Gray)

The color in parentheses is the color of the cartridge attached to the stapler. It is not the color of NEOVEIL™.

NEOVEIL™ Sheet
Product video



NEOVEIL™ SLR
Product video



Distributor

GUNZE
GUNZE MEDICAL

GUNZE LIMITED Medical Division

<http://www.gunze.co.jp/e/medical/>

46, Natsumegaichi, Aono-cho, Ayabe,
Kyoto 623-8513, Japan

U.S.A.

GUNZE INTERNATIONAL USA, INC.

390 Fifth Avenue, Suite 903, New York,
NY 10018, U.S.A.

Tel: +1-212-354-9060

Fax: +1-212-354-6171

inquiry_gius@gunze.co.jp

<https://gunzemedicalusa.com/>

Europe

**GUNZE INTERNATIONAL
EUROPE GmbH**

Oststraße 80, 40210 Düsseldorf, Germany

Tel: +49-211-3613-971

Fax: +49-211-3613-164

China

**GUNZE MEDICAL DEVICES
(SHENZHEN) LIMITED**

Room910-911, Material Landmark Building,
NO.6 Liyuan Road, Sungang Street,
Luohu District Shenzhen 518022, P.R. China

Tel: +86-755-8230-0553

Fax: +86-755-8230-0469

Japan/Asia Pacific

GUNZE MEDICAL LIMITED TOKYO OFFICE

15F, Tokyo Shiodome Bldg., 1-9-1,
Higashi Shimbashi, Minato-ku,
Tokyo 105-7315, Japan

Tel: +81-3-4485-0020

Fax: +81-3-4485-0050